
Calyx RIM 7.1.1 Release Notes

CALYX™

Contents

Calyx RIM 7.1.1 Release Notes.....	1
Calyx RIM 7.1.1 Release Notes.....	1
Calyx RIM 7.1.1 Resolved Issues.....	1
Calyx RIM7.1 COB and Known Issues.....	2
Index.....	3

Calyx RIM 7.1.1 Release Notes

Calyx RIM 7.1.1 Release Notes

This release is not a general release of Calyx RIM 7.1.1.

This release includes the following database scripts.

Component	Customizable	Build Number
Migration	No	DB_update7_1_0_0000_0114
Post-Migration	No	DB_update7_1_1_0000_0001

Note: The `DB_update7_1_0_0000_0114` migration script resolves all the errors listed in the *Calyx RIM 7.1.1 Resolved Issues* topic. This migration script must be used when migrating. The *Post-Migration* script resolves only the 26961 defect listed in the *Calyx RIM 7.1.1 Resolved Issues* topic and must be used if you have already migrated to 7.1.

Calyx RIM 7.1.1 Resolved Issues

The following 7.1 migration issues have been reported and have been resolved.

Issue	Description
27305	ORA-54033 errors occurred when using the <code>DB_update7_1_0_0000_0112</code> database migration script to migrate to 7.1 if some of the <code>CREATION_DATE/BEG_DATE</code> fields were defined to use the <code>DATE</code> data type instead of the <code>TIMESTAMP</code> data type. This issue is resolved when migrating to 7.1 using the <code>DB_update7_1_0_0000_0114</code> database migration script included in this release.
27304	Having a product with a very large number of applications and full product presentations resulted in an <code>ORA-01489</code> error when using the <code>DB_update7_1_0_0000_0112</code> database migration script to migrate to 7.1. This issue is resolved when migrating to 7.1 using the <code>DB_update7_1_0_0000_0114</code> database migration script included in this release.
26961	Substances in the PDS tree were not displayed when using the <code>DB_update7_1_0_0000_0112</code> database migration script to migrate to 7.1. This issue is resolved when migrating to 7.1 using the <code>DB_update7_1_0_0000_0114</code> database migration script included in this release. For customers that have already migrated successfully to 7.1, this issue can be resolved by using the <code>DB_update7_1_1_0000_0001</code> post-migration script included in this release.

Calyx RIM7.1 COB and Known Issues

For Clarification of Behavior (COB) and Known Issues, see: COB and Known Issues topics under Calyx RIM7.1 Release Notes (Cumulative).

Index

R

Release Notes [2](#)

Resolved Issues [1](#)