
Calyx RIM 7.1.13 Release Notes

CALYX™

Contents

Calyx RIM 7.1.13 Release Notes	1
Calyx RIM 7.1.13 Resolved Issues.....	1
Calyx RIM7.1 COB and Known Issues.....	1
Index	2

Calyx RIM 7.1.13 Release Notes

The Release Notes contain information about updated files, resolved issues, and any open issues or workarounds in this release of Calyx RIM 7.1.13.

The details in these Release Notes supersede all similar details that appear anywhere else in this release of the product or its other documentation.

Calyx RIM 7.1.13 Resolved Issues

Calyx RIM 7.1.13 resolves the following issues.

Publishing

Issue	Description
DE31334, [TASK0444686]	When Recalculate Checksum was performed on an assembly with 1000+ leafs, the Job Requests > Job Message query displayed the <code>Error: could not execute query</code> message. This issue is resolved. Now, the Recalculate Checksum operation completes correctly and the Job Message displays relevant information.
DE31296, [TASK0373863]	When working with large assemblies, performing the Recalculate Checksum operation even on a single element timed out and failed. This issue is resolved. Now, the performance is improved, and the Recalculate Checksum operation is performed only on the files selected by the user.

Calyx RIM7.1 COB and Known Issues

For Clarification of Behavior (COB) and Known Issues, see: [COB](#) and [Known Issues](#) topics under Calyx RIM7.1 Release Notes (Cumulative).

Index

N

New Features [1](#)

R

Release Notes [1](#)

Resolved Issues [1](#)