
Calyx RIM 7.1.2 Release Notes

CALYX™

Contents

Calyx RIM 7.1.2 Release Notes.....	1
Calyx RIM 7.1.2 Release Notes.....	1
Calyx RIM 7.1.2 Resolved Issues.....	1
Calyx RIM7.1 COB and Known Issues.....	1
Index.....	2

Calyx RIM 7.1.2 Release Notes

Calyx RIM 7.1.2 Release Notes

This document describes a sub-minor release provided to fix the defects found with Calyx RIM version 7.1. The Quality Control team performed first product install and regression testing on this sub-minor release.

Calyx RIM 7.1.2 Resolved Issues

Calyx RIM 7.1.2 resolves the following issues.

Issue	Build	Description
DE27634; [INC0084824]	7.1.2.0.0013	When running the <i>Global Project Plan</i> wizard with the Create CTA/IND Applications and Events project type selected and clicking Next on <i>Select the protocol</i> , the <i>Enter Application properties</i> page timed out. The wizard could not load the Applications. This performance issue is resolved. The <i>Global Project Plan</i> wizard loads Applications successfully.
DE27552; [INC0084364]	7.1.2.0.0013	Publishing a Documentum document with a null keyword attribute value resulted in the publish job failing with the error <code>Unable to get DMS attribute 'keywords[]' for document XXX</code> . This issue is resolved.
DE27527; [INC0082927, INC0085354]	7.1.2.0.0013	When publishing STFs using an Output Folder setting along with Create Study Folders = No , the published output contained 2 STF Folders, one at the root containing the STF.xml and another in the correct location containing the associated PDF. This issue is resolved.
DE27515; [INC0082577]	7.1.2.0.0013	Performance issues were experienced in the following cases: <ul style="list-style-type: none"> – After the Withdrawal Workflow to retrieve Registrations after Application entered. – When Registration was created from the Registrations tab on the Application page. – During FPP loading. These issues are resolved. The performance of the queries is now improved.

Calyx RIM7.1 COB and Known Issues

For Clarification of Behavior (COB) and Known Issues, see: COB and Known Issues topics under Calyx RIM7.1 Release Notes (Cumulative).

Index

R

Release Notes 1
Resolved Issues 1