

---

# Calyx RIM 7.1.3 Release Notes

# Contents

<b>Calyx RIM 7.1.3 Release Notes.....</b>	<b>1</b>
Calyx RIM 7.1.3 Release Notes.....	1
Entity Fields.....	1
Fields and Columns Added: Calyx RIM 7.1.3.....	1
Database Columns Changed: Calyx RIM 7.1.3 .....	2
Resolved Issues: Calyx RIM 7.1.3.....	3
Calyx RIM7.1 COB and Known Issues.....	4
<b>Index.....</b>	<b>5</b>

## Calyx RIM 7.1.3 Release Notes

### Calyx RIM 7.1.3 Release Notes

The Release Notes contain information about updated files, new functionality, resolved issues, and any open issues or workarounds in this release of Calyx RIM 7.1.3.

The details in these Release Notes supersede all similar details that appear anywhere else in this release of the product or its other documentation.

#### Regulatory Updates

Support for China (CN) eCTD Module 1 version 1.0.

#### API Updates

Calyx RIM 7.1.3 includes multiple API updates including updates for IDMP.

API support for Application Country EU CTR Applications. Any entity created or updated as a result of the EU CTR support is API enabled.

Changes between Calyx RIM 7.1 API and Calyx RIM 7.1.3 API are documented in the *Calyx RIM Administration*. See: *Calyx RIM Administration > Technical Administration > Calyx RIM API Documentation*.

#### IDMP

Additional IDMP data points for capturing alternate active ingredients for Pharmaceutical Products and Manufactured Items.

#### Power BI

Support for the SQL Direct Access model that includes the SQL reporting database with support for the Registration, Security Data and tables that support out of the box Power BI reports.

### Entity Fields

#### Fields and Columns Added: Calyx RIM 7.1.3

\

#### Entities

The following table describe fields and columns that were added in this release. These include the fields and column that were added to support IDMP data capture.

Table 1: Active Ingredients

Field/Option Name	Field Type	Description	Business Rules
Alternate IDMP Active Ingredient	Drop-down list, selection widget.	The <i>Pharmaceutical Product Active Ingredient Attributes</i> and <i>Manufactured Item Active Ingredient Attributes</i> pages now include the respective field. See: <i>Registrations &gt; Active Ingredient Attributes</i> .	<ul style="list-style-type: none"> <li>For Manufactured Item Active Ingredient: this field can be configured as a drop-down list with the values from the Component Reference Active Ingredients.</li> <li>For Pharmaceutical Product Active Ingredient: this field can be configured as a drop-down list with the values from the Component Reference Active Ingredients of Components associated with a Pharmaceutical Product. To search for a Component Active Ingredient, you need to assign Components to Pharmaceutical Product.</li> </ul> <p>This field can also be configured as the substance selection widget with all the active values from <b>Data Administration &gt; Substance Values</b>.</p> <hr/> <p><i>Note: For configuration assistance, please contact your Business Development Representative.</i></p> <hr/>

### Database Columns Changed: Calyx RIM 7.1.3

The following Database Columns were added, modified, or removed in this release.

Field Name	Field Size	Database Column Name	Db Column Size	Db Table Name
Alternate IDMP Active Ingredient	N/A	ALTERNATE_SUBSTANCE_ID	NUMBER	XPHARMA_PRODUCT_AI
	N/A	ALTERNATE_SUBSTANCE_ID	NUMBER	MI_AI

## Resolved Issues: Calyx RIM 7.1.3

The following issues that were reported to Technical Support have been resolved. Issues are listed in descending order by their tracking numbers. Customer incident numbers that were recorded when the issues were reported are included in brackets [#####] following the tracking number.

If an issue was resolved in a previously released certified hotfix (CHF), that CHF build number is included.

### Publishing

Issue	Build	Hotfix Build	Description
DE27672 [INC0085185]	7.1.3.0.0006	N/A	When searching for a document based on EDMS Number in <i>DMS Browse</i> window to assign it to the assembly, it was not displayed and no custom attributes were populated. This issue is resolved. The document with the custom attributes is displayed.
DE27588; [INC0084413]	7.1.3.0.0006	N/A	When an assembly contained docrefs that were in a repository that was not specified in <i>Technical Administration</i> , the import of an assembly failed. This issue is resolved. The assembly imports with the documents unassigned.
DE26892	7.1.3.0.0006	N/A	China: After importing an eCTD folder, the 2.3.1 section was missing in the structure of the imported assembly. This was caused by incorrect published output. This issue is resolved. Now, the structure of the imported assembly will be the same as the source. The CN eCTD ICH Module2-5 v3.2.xml file template is updated for China region with this fix.
26825	7.1.3.0.0006	6.2.0.8	EAEU: If the Drug Registration Document Code field had a value in a submission, the Drug Registration Document Name was not included in the index-r-022.xml file after publishing. The Drug Registration Document Name was also not saved in the leaf after an eCTD import. The issue is resolved. Now, both Drug Registration Document Name and Drug Registration Document Code are included in the index-r-022.xml file after publishing and in the leaf after an eCTD import.

**Registrations**

Issue	Build	Hotfix Build	Description
26768; [INC0058320]	7.1.3.0.0006	6.2.0.8	<p>The /all and /query endpoints did not work via API for the package set registration, product registration, full product presentation, package set - country entities and resulted in an error. This issue is resolved. The following endpoints work correctly via API:</p> <ul style="list-style-type: none"> <li>– /insight/api/v2/license/packageset-registration/all</li> <li>– /insight/api/v2/license/product-registration/all</li> <li>– /insight/api/v2/license/full-product-presentation/all</li> <li>– /insight/api/v2/license/packageset/all</li> <li>– /insight/api/v2/license/packageset/country/all</li> </ul>

**Calyx RIM7.1 COB and Known Issues**

For Clarification of Behavior (COB) and Known Issues, see: COB and Known Issues topics under Calyx RIM7.1 Release Notes (Cumulative).

# Index

## N

New Features [1, 2](#)

## R

Release Notes [4](#)

Resolved Issues [3](#)