
Calyx RIM 7.1.6 Release Notes

CALYX™

Contents

Calyx RIM 7.1.6 Release Notes.....	1
Calyx RIM 7.1.6 Release Notes.....	1
Calyx RIM 7.1.6 Resolved Issues.....	1
Calyx RIM7.1 COB and Known Issues.....	1
Index.....	2

Calyx RIM 7.1.6 Release Notes

Calyx RIM 7.1.6 Release Notes

The Release Notes contain information about updated files, new functionality, resolved issues, and any open issues or workarounds in this release of Calyx RIM 7.1.6.

The details in these Release Notes supersede all similar details that appear anywhere else in this release of the product or its other documentation.

This release resolves issues with missing database indexes.

Note: Calyx RIM 7.1.6 is a database-only release and includes the DB_update7_1_6_0000_0001.zip file only.

Calyx RIM 7.1.6 Resolved Issues

Calyx RIM 7.1.6 resolves the following issue.

Publishing

Issue	Build Number	Description
DE28488; [INC0108112]	7.1.5.0.0004 + DB_update7_1_6_0000_0001	Multiple users creating/updating different entities simultaneously resulted in Calyx RIM being unusable for a period of time. This issue is resolved.

Calyx RIM7.1 COB and Known Issues

For Clarification of Behavior (COB) and Known Issues, see: COB and Known Issues topics under Calyx RIM7.1 Release Notes (Cumulative).

Index

N

New Features [1](#)

R

Release Notes [1](#)

Resolved Issues [1](#)