

---

# Calyx RIM 7.1.9 Release Notes

CALYX™

## Contents

|   |          |
|---|----------|
| <b>Calyx RIM 7.1.9 Release Notes.....</b>   | <b>1</b> |
| Calyx RIM 7.1.9 Release Notes.....          | 1        |
| Entity Fields Updated: Calyx RIM 7.1.9..... | 1        |
| Calyx RIM 7.1.9 Resolved Issues.....        | 2        |
| Calyx RIM7.1 COB and Known Issues.....      | 2        |
| <b>Index.....</b>                           | <b>3</b> |

## Calyx RIM 7.1.9 Release Notes

### Calyx RIM 7.1.9 Release Notes

The Release Notes contain information about updated files, resolved issues, and any open issues or workarounds in this release of Calyx RIM 7.1.9.

The details in these Release Notes supersede all similar details that appear anywhere else in this release of the product or its other documentation.

#### Regulatory Updates

EAEU v1.0:

- The index-r-022.xml file now includes new checksum attributes to conform to EAEU v1.0 eCTD requirements.
- Publishing functionality is updated with the checksum information in index-r-022.xml for each document. The checksum attributes are removed for EAEU XML with embedded PDFs during the export and import of EAEU submissions.
- The Export Publish Output page is updated to get the copy of submission without checksums with the help of the XML without checksum option.

#### Entity Fields Updated: Calyx RIM 7.1.9

Entities fields updated in Calyx RIM 7.1.9.

The following entity field is updated in Calyx RIM 7.1.9.

| Field Name                | Field Type | Entity Name  | Description | Business Rules  |
|---------------------------|------------|--|-------------|---|
| Embed PDFs into EAEU XML? | Box        | This field is renamed to EAEU Output and updated to the list with the Embed PDFs and XML without checksum options on the Export Published Output page. | EAEU Output | You can use EAEU Output to export EAEU submissions. See: <i>Calyx RIM Best Practice: EAEU Submissions</i> . |

## Calyx RIM 7.1.9 Resolved Issues

Calyx RIM 7.1.9 resolves the following issues.

### Publishing

| Issue                   | Description  |
|-------------------------|--|
| DE29280<br>[INC0145064] | <p>Adding the "/" (forward slash) character in assembly metadata resulted in the following issues after publishing.</p> <ul style="list-style-type: none"> <li>- Empty attributes in published output XML for 2.3.P Leaf.</li> <li>- Missing galenic-name XML elements in the ch-regional.xml file.</li> </ul> <p>Both issues are now resolved. Now, assembly metadata is published correctly even with the "/" (forward slash) character added.</p> |
| DE29254                 | <p>In an assembly with empty leafs, running the Recalculate checksum operation failed. This issue is resolved. Now, the Recalculate checksum operation completes.</p>  |
| DE29174<br>[INC0143854] | <p>When importing an assembly with assigned source document from assembly file, the source location path in the <b>Document Attributes</b> window was not populated. This issue is resolved. Now, the source location displays the location path to the assigned source document.</p>  |
| DE28851<br>[INC0123232] | <p><b>Group Submissions:</b> In an <b>Assembly Plan</b> with multiple applications, when two folders containing leafs and a form were added for each application, the xml file displayed leaf elements under one application entry after publishing. This issue is resolved. Now, the xml file displays leaf elements under each application entry.</p>  |

## Calyx RIM7.1 COB and Known Issues

For Clarification of Behavior (COB) and Known Issues, see: COB and Known Issues topics under Calyx RIM7.1 Release Notes (Cumulative).

# Index

## N

New Features [1](#)

## R

Release Notes [1](#), [2](#)

Resolved Issues [2](#)