

e|NNOV

Contents

Ennov InSight 7.2.10 Release Notes	1
Ennov InSight 7.2.10 Release Notes.....	1
Resolved Issues: Ennov InSight 7.2.10	1
Ennov InSight 7.2 COB and Known Issues.....	2
Data Administration Tables	3
Data Administration Sections Added: Ennov InSight 7.2 .10.....	3
Index	4

Ennov InSight 7.2.10 Release Notes

Ennov InSight 7.2.10 Release Notes

The Release Notes contain information about updated files, new functionality, resolved issues and any open issues or workarounds in this release of Ennov InSight 7.2.10.

The details in these Release Notes supersede all similar details that appear anywhere else in this release of the product or its other documentation.

Regulatory Updates

Support for South Africa (ZA) eCTD Module 1 version 3.0.

Note: This functionality is for presentation purpose only. ZA 3.0 is not supported by the agency and therefore cannot be used for a submission as it won't be accepted.

For more information about new ZA values, extensions and rules, see:

- *Ennov InSight 7.2 Publisher > Submissions > Create eCTD Wizard > eCTD Attributes South Africa*
- *Ennov InSight 7.2 Publisher > Extension Type and Usage > Extension Type: Regional Folders > South Africa*
- *Ennov InSight 7.2 Publisher > Administrative Activity > Import eCTD Applications > Prepare for Bulk eCTD Import > ZA 3.0 Applications*
- *Ennov InSight 7.2 Administration > eCTD Codes > eCTD Filing Types for ZA*
- *Ennov InSight 7.2 Administration > eCTD Filing Types for ZA > eCTD Sub-Filing Types for ZA*
- *Ennov InSight 7.2 Administration > Assembly XML Specification > m1-za v. 3.0*

Resolved Issues: Ennov InSight 7.2.10

The following issues that were reported to Technical Support have been resolved. Issues are listed in descending order by their tracking numbers. Customer incident numbers that were recorded when the issues were reported are included in brackets [#####] following the tracking number.

If an issue was resolved in a previously released sub- minor build, that build number is included.

Common

Issue	Sub-Minor Build	Description
DE31815, [TASK0925392]	N/A	If Documentum Virtual Document (VD) with child elements was assigned to an assembly, the display order of assembly elements was not in the same order as in Documentum. This issue is resolved. Now, the assigned VD retains the original order defined in Documentum.
DE31799, [INC0321912, TASK0869957]	N/A	Trying to assign a virtual document using DMS Search resulted in a placeholder being displayed instead of the expected content. This issue is resolved. Now, a virtual document assigned using DMS Search is displayed correctly.

Addendum 1. Certified on Ennov InSight 7.2.8

Issue	Sub-Minor Build	Description
DE31757; [INC0321912, TASK0869957]	N/A	If the Documentum file names included invalid characters: <code>\\/:*?\"<> </code> , performing some actions in assembly resulted in an error. This issue is resolved. Invalid symbols in file names are now automatically deleted, and there are no errors when working with assembly.

Addendum 2. Certified on Ennov InSight 7.2.8

Issue	Sub-Minor Build	Description
DE31795	N/A	Missing indexes in the database caused performance issue when trying to save entities after editing. This issue is resolved. Run <code>DB_update7_2_(7-8-9)_DE31795_0001.zip</code> script in addition to the main scripts during the migration procedure. The additional script must be run after the main script. See <i>Ennov InSight Database Migration Details</i> .

Ennov InSight 7.2 COB and Known Issues

For Clarification of Behavior (COB) and Known Issues, see: [COB](#) and [Known Issues](#) topics under Ennov InSight 7.2 Release Notes (Cumulative)

Data Administration Tables

Data Administration Sections Added: Ennov InSight 7.2 .10

The following Data Administration tables, sections, options, widgets, filters or fields were added for Ennov InSight 7.2.10.

Data Administration Section: Submission Maintenance

Table 1: Evaluation Pathway Values

Entity Affected	Field Name	Value List Description
Assembly	<ul style="list-style-type: none"> – Region – Applicable DTD/Schema – eCTD Code – Evaluation Pathway – Last User Updated – Last Changed Date – Active Flag 	To conform to the latest standards set by South African Health Products Regulatory Authority (SAHPRA), the <i>Evaluation Pathway Values</i> table is introduced.

Index

N

New Features [1](#), [3](#)

R

Release Notes [2](#)

Resolved Issues [1](#)