



Contents

Ennov InSight 7.2.11 Release Notes	1
Ennov InSight 7.2.11 Release Notes.....	1
Resolved Issues: Ennov InSight 7.2.11	1
Ennov InSight 7.2 COB and Known Issues.....	2
Index	3

Ennov InSight 7.2.11 Release Notes

Ennov InSight 7.2.11 Release Notes

The Release Notes contain information about updated files, new functionality, resolved issues and any open issues or workarounds in this release of Ennov InSight 7.2.11.

The details in these Release Notes supersede all similar details that appear anywhere else in this release of the product or its other documentation.

Support of Special Characters

Ennov InSight 7.2.11 supports the use of more special characters in Documentum and File system. *See: Ennov InSight Administration 7.2 > Special Characters.*

Regulatory Updates

Support for South Africa (ZA) eCTD Module 1 version 3.1.

For more information about new ZA values, extensions, and rules, see:

- *Ennov InSight 7.2 Data Migration > Ennov InSight 7.2.11 Updates*
- *Ennov InSight 7.2 Publisher > Submissions > Create eCTD Wizard > eCTD Attributes South Africa*
- *Ennov InSight 7.2 Publisher > Extension Type and Usage > Extension Type: Regional Folders > South Africa*
- *Ennov InSight 7.2 Publisher > Administrative Activity > Import eCTD Applications > Prepare for Bulk eCTD Import > ZA 3.1 Applications*
- *Ennov InSight 7.2 Administration > eCTD Codes > eCTD Filing Types for ZA*
- *Ennov InSight 7.2 Administration > eCTD Filing Types for ZA > eCTD Sub-Filing Types for ZA*
- *Ennov InSight 7.2 Administration > Assembly XML Specification > m1-za v. 3.1*

Note: *As the South Africa (ZA) eCTD Module 1 version 3.0 is not accepted by SAHPRA, it is no longer supported in Ennov InSight .*

Resolved Issues: Ennov InSight 7.2.11

The following issues that were reported to Technical Support have been resolved. Issues are listed in descending order by their tracking numbers. Customer incident numbers that were recorded when the issues were reported are included in brackets [#####] following the tracking number.

If an issue was resolved in a previously released sub- minor build, that build number is included.

Common

Issue	Sub-Minor Build	Description
DE31877	N/A	For a Documentum Virtual Document (VD) assigned to an assembly, changes made to its structure (for example, adding a new child document) were not displayed after a VD was reassigned. This issue is resolved. Now, a VD structure correctly displays the changes after a VD is reassigned.

Ennov InSight 7.2 COB and Known Issues

For Clarification of Behavior (COB) and Known Issues, see: [COB](#) and [Known Issues](#) topics under Ennov InSight 7.2 Release Notes (Cumulative)

Index

N

New Features [1](#)

R

Release Notes [2](#)

Resolved Issues [1](#)